



## Configurations Delta 1-24 / 2-24

- 1 freezing and drying inside ice condenser chamber (process A) on 10 temperature controlled shelves ( $\varnothing$  200 mm, total area = 0.31 m<sup>2</sup>, distance between shelves = 25 mm - can be varied by removing single shelves)
- 2 freezing and drying inside ice condenser chamber (process A) on 4 temperature controlled shelves with stoppering device ( $\varnothing$  250 mm, total area = 0.18 m<sup>2</sup>, distance between shelves = 45 mm - can be increased to 110 mm by removing one shelf)
- 3 2 drying chambers with 24 ports for flasks, wide-neck filter bottles or distributors for ampoules
- 4 10 heatable shelves ( $\varnothing$  200 mm, total area = 0.31 m<sup>2</sup>, distance between shelves = 25 mm - can be varied by removing single shelves) additional ports for 12 round bottom flasks, wide-neck filter bottles or distributors for ampoules
- 5 4 heatable shelves with stoppering device ( $\varnothing$  250 mm, total area = 0.18 m<sup>2</sup>, distance between shelves = 50 mm - can be varied by removing single shelves) additional ports for 12 round bottom flasks, wide-neck filter bottles or distributors for ampoules
- 6 8 heatable shelves ( $\varnothing$  375 mm, total area = 0.88 m<sup>2</sup>, distance between shelves = 48 mm - larger distances can be obtained by special order)



T 0346 26 90 90  
F 0346 26 48 14

I www.salmenkipp.nl  
E info@salm-en-kipp.nl