

Configurations Epsilon 1-4 / 2-4

A variety of Epsilon 1-4 / 2-4 configurations are possible due to a very clever and variable system of shelves that serve as an intermediate valve and defrosting plate for the ice condenser. A selection of these combinations is illustrated here.

- 1 base system with combination stainless steel dish for product or for defrosting-water
- 2 system with stoppering device and second shelf (intermediate valve) plus integrated manifold
- 3 illustration of closing vials



salmenkipp

T 0346 26 90 90
F 0346 26 48 14

I www.salmenkipp.nl
E info@salm-en-kipp.nl