

Pipetman® F

Fixed-Volume Pipettes

Pipetman F pipettes are fixed-volume air-displacement pipettes. Thirteen models cover the volume range from 2 µL to 1000 µL.

◎ The right pipette

- No risk of volume-selection errors. Using a dedicated fixed-volume pipette for a specific test assures reliable and consistent results.

◎ The right volume

- Volume can be adjusted by users in the laboratory for dense or viscous fluids (once reset, the new volume cannot be changed accidentally).

◎ The right price

- The single-volume mechanism provides a low-cost alternative for clinical diagnostics, quality control and routine testing.

Pipetman F Models				
Model	Diamond Tips	Part Number	Volume (µL)	List Price
F2	D200	F123770	2	94,80 €
F5	D200	F123771	5	94,80 €
F10	D200	F123772	10	94,80 €
F20	D200	F123604	20	94,80 €
F25	D200	F123775	25	94,80 €
F50	D200	F123778	50	94,80 €
F100	D200	F123784	100	94,80 €
F200	D200	F123605	200	94,80 €
F250	D1000	F123787	250	94,80 €
F300	D1000	F123788	300	94,80 €
F400	D1000	F123789	400	94,80 €
F500	D1000	F123790	500	94,80 €
F1000	D1000	F123606	1000	94,80 €

Gilson maximum permissible errors are guaranteed only when Pipetman F pipettes are used with the recommended Gilson Diamond tips.

Distriman®

Repetitive Pipette

Distriman is a positive-displacement, continuously adjustable repeater-pipette that simplifies multiple dispensing when used with Gilson pure polypropylene DistriTip® syringes.

◎ Direct volume readout

- Dispense volume reads directly in microliters, eliminating troublesome calculations.



Direct volume readout

◎ Only three syringes are needed

- Micro 125 µL, Mini 1250 µL and Maxi 12.5 mL syringes cover the entire volume range (from 1 µL to 1.25 mL). DistriTips syringes are the only syringes capable of dispensing up to 100 aliquots.

◎ Selection of packaging

- Each DistriTips syringe model is available in standard or sterilized (ST) form, in boxes of 50. Sterilized DistriTips are individually wrapped and certified free of RNase, DNase, DNA, RNA and protease.



Selection of packaging

Distriman Repetitive Pipette				
	Part Number		Price	
	F164001		219,60 €	
DistriTips Syringes (box of 50)				
Model	Aliquot Range	Part Number	Volume (µL)	Price
125 µL	1 µL to 12.5 µL	F164100	Min.	2
			Max.	10
Micro Micro ST†		F164130	Min.	5 33,00 €
			Max.	10 67,98 €
1250 µL	10 µL to 125 µL	F164110	Min.	20
			Max.	50 33,00 €
Mini Mini ST		F164140	Min.	200
			Max.	100 67,98 €
12.5 mL	100 µL to 1.25 mL	F164120	Min.	200
			Max.	500 33,00 €
Maxi Maxi ST		F164150	Max.	1000 67,98 €

†ST means sterilized syringes (individually wrapped).

