

# Precise Dosing and Dispensing

Premium Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Safety for continuous operation: the motor will be switched off if a high thermal load situation occurs to **prevent accidents**

**High precision** even for flow rates starting at 0.005 ml/min



Optional foot-pedal remote control via cable to operate pump drive **outside a closed fume hood** facilitates operation

All models feature a smooth start operation which **prevents spills and splashing media**. The speed will ramp up slowly until your set rpm has been reached

For the use in cell biology just use a pump head which features convex rollers and **does not squeeze** the tubing as conventional rollers

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a truly worthwhile investment.

**salmenkipp**

T 0346 26 90 90  
F 0346 26 48 14

I www.salmenkipp.nl  
E info@salm-en-kipp.nl