

PIPETMAN® Neo

The Legend is Reborn!



Your benefits

PREMIUM ROBUSTNESS

- Long product life and high-resistance to chemical attack with its stainless-steel and PVDF materials
- Air-tight seal with its unique "dry seal" technology (grease-free)

ACCURACY & PRECISION

- Performance exceeding ISO recommendations by an average of 50%

COMFORTABLE OPERATION

- Increased pipetting comfort with lowered pipetting forces

Gilson Secured Pipetting System

ISO 8655-2 recommends the use of tips supplied by your pipette manufacturer for optimal results. DIAMOND tips guarantee that PIPETMAN Neo meets our specified tolerances.

Complete your equipment with



Trio stand



The Jimmy



0,2µL - 1000µL



Non viscous
Non volatile
Non dense

salmenkipp

T 0346 26 90 90
F 0346 26 48 14

I www.salmenkipp.nl
E info@salm-en-kipp.nl