

Vacuum Safety Suction Systems - the correct selection



The Automaticals

Art.-No.	Article
AA 02	Vacuum Safety Suction System, with 2 litre bottle, complete with pump, tube, filter etc.
AA 04	Vacuum Safety Suction System, with 4 litre bottle, complete with pump, tube, filter etc.

With automatic vacuum pump, suction capacity: 30 l / min (air)



The Compacts

Art.-No.	Article
AC 02	Vacuum Safety Suction System, with 2 litre bottle, complete with pump, tube, filter etc.
AC 04	Vacuum Safety Suction System, with 4 litre bottle, complete with pump, tube, filter etc.

With a powerful membrane pump, suction capacity: 15 l / min (air)



The Simples

Art.-No.	Article
AZ 02	Vacuum Safety Suction System, with 2 litre bottle, complete with pump, tube, filter etc.
AZ 04	Vacuum Safety Suction System, with 4 litre bottle, complete with pump, tube, filter etc.

With a small, silent membrane pump, suction capacity: 4 l / min (air)

Pumps / Accessories for Suction Systems

Pumps

Art.-No.	Article
AP 04	Membrane Vacuum Pump, 4 l / min (air), 250 mbar
AP 15	Membrane Vacuum Pump, 15 l / min (air), 100 mbar
AF 01	Foot-operated switch

Comprehensive parts and options

AE 02	Safety Collection Unit, with 2 litre bottle, complete, but without pump
AE 04	Safety Collection Unit, with 4 litre bottle, complete, but without pump
AL 02	Spare bottle 2 litre, with fast locks in screw cap
AL 04	Spare bottle 4 litre, with fast locks in screw cap
AW 01	Hose clip, in addition to the spare bottles
AS 02	Safety Bottle Holder for 2 litre bottle, anodised aluminium
AS 04	Safety Bottle Holder for 4 litre bottle, anodised aluminium
AV 03	Safety Filter, packing unit: 3 pcs.

Handles and adapters

AM 01	Handle for quantity regulation, for absorption of a Pasteur pipette or similar
AH 01	Handle with trumpet valve, with adapter for absorption of a Pasteur pipette or similar
AR 05	Stainless steel suction adapter, 50 mm long, special for AH 01
AR 15	Stainless steel suction adapter, 150 mm long, special for AH 01
AR 25	Stainless steel suction adapter, 250 mm long, special for AH 01
AR 40	Suction adapter for yellow pipette tips, special for AH 01
AR 48	Adapter with suction rake, 8-canal, for titer plates, autoclavable, special for AH 01
AR 58	Suction rake, 8-canal, for titer plates, autoclavable, as spare part
AY 01	Y-Distributor for connecting several collection units

salmenkipp

T 0346 26 90 90 I www.salmenkipp.nl
F 0346 26 48 14 E info@salmenkipp.nl