

REMEMBER YOUR QUOTE ON OPTIONALS

LTF / LCF series

As a synonym for security, the LTF and LCF series gives reliable low temperature freezing, with good temperature uni-formity. Digital display, alarms, contact for remote monitoring, basket(s) and castors are some of the standards.

LTF / LCF series

- Single compressor
- Low energy consumption
- Low noise level
- 100% HCFC/CFC free
- Key lock, castors and baskets



LTF controller



Controller features

- Digital display
- Visual and acoustic alarm
- Adjustable high/low temperature alarm
- Contact for remote alarm
- Probe failure alarm



- Remember your quote on:
- Racks and boxes
 - Chart Recorder
 - CO2 Backup
 - GSM Alarm

Model	LTF 85	LTF 325	LTF 425	LTF 535
Item code (DAI no.)	DAI101	DAI111	DAI116	DAI121
Dimensions exterior WxDxH (mm)	552x648x850	1262x698x885	1562x698x885	1662x758x925
Net capacity (L)	71	284	368	504
Temperature range (°C)	-30/-60	-30/-60	-30/-60	-30/-60
Power supply (V)	230	230	230	230
Frequency (Hz)	50/60	50/60	50/60	50/60

Model	LCF 100	LCF 200	LCF 300	LCF 400	LCF 500
Item code (DAI no.)	DAI177	DAI178	DAI179	DAI180	DAI181
Dimensions exterior WxDxH (mm)	725x730x865	1055x730x865	1305x730x865	1505x730x865	1705x730x865
Net capacity (L)	130	230	300	360	485
Temperature range (°C)	-30/-60	-30/-60	-30/-60	-30/-60	-30/-60
Power supply (V)	230	230	230	230	230
Frequency (Hz)	50	50	50	50	50



We reserve the right to change specifications without notice. Subject to confirmation and availability.