

# Microman®

## The World Standard for Positive-Displacement Pipettes

Gilson's Microman® is a fully-adjustable, positive-displacement pipette that, together with its disposable Capillary Pistons (CPs), will assure accuracy and precision—even with problem liquids.

### ⊙ Improve accuracy and reliability with problem liquids

- The perfect solution for complete and rapid pipetting of viscous and dense liquids such as oil, syrup, cosmetic cream, liquified food, paint, glycerol or buffers.
- Accurately measure high-vapor-pressure liquids such as acetone, chloroform, alcohol, hexane and other solvents without leakage or dribbling.

### ⊙ 100% protection against contamination

- Disposable capillaries and pistons isolate the pipetting instrument from samples and aerosols, thereby preventing sample cross-contamination.
- Built-in push-button tip-ejector prevents user contact with contaminated capillaries and pistons.
- Contamination-free pipetting with Microman pipettes helps prevent false results in PCR analysis.

### ⊙ Microman and Capillary Pistons (CPs)

- Volume range from 1 µL to 1000 µL with six different models.
- Color-coded packaging identifies correct pipette model.
- No need for recalibration when you change the CP.
- Disposable, safe – Capillaries made of pure, unbreakable polypropylene.
- Dye-free – no risk of chemical contamination.
- Slender profile – Microman CPs fit even the narrowest and deepest tubes.

Microman Pipettes				
Model	Capillary Piston <sup>†</sup>	Part Number	Volume (µL)	Price
M10 1 - 10 µL	CP10	F148501	Min. 1 5 Max. 10	185,25 €
M25 3 - 25 µL	CP25	F148502	Min. 3 10 Max. 25	185,25 €
M50 20 - 50 µL	CP50	F148503	Min. 20 Max. 50	185,25 €
M100 10 - 100 µL	CP100	F148504	Min. 10 50 Max. 100	185,25 €
M250 50 - 250 µL	CP250	F148505	Min. 50 100 Max. 250	185,25 €
M1000 100 - 1000 µL	CP1000	F148506	Min. 100 500 Max. 1000	185,25 €



Microman Capillary Pistons							
Capillary Piston <sup>†</sup>	Microman Model	Part Number	Quantity/Box	Sterilized	Assembled	Racked in Tipack	Price
CP10	M10	F148412	192	No	Yes	Yes	52,90 €
CP10	M10	F148312	960	No	Yes	Yes	221,75 €
CP25	M25	F148112	200	No	No	No	49,15 €
CP50	M50	F148113	200	No	No	No	49,15 €
CP100	M100	F148414	192	No	Yes	Yes	52,90 €
CP100	M100	F148314	960	No	Yes	Yes	221,75 €
CP250	M250	F148114	200	No	No	No	49,15 €
CP1000	M1000	F148560	182	No	Yes	Yes	52,90 €
CP10ST	M10	F148413	192	Yes	Yes	Yes	59,85 €
CP10ST	M10	F148313	960	Yes	Yes	Yes	252,00€
CP100ST	M100	F148415	192	Yes	Yes	Yes	59,85 €
CP100ST	M100	F148315	960	Yes	Yes	Yes	252,00€
CP1000ST	M1000	F148180	182	Yes	Yes	Yes	59,85 €

<sup>†</sup>Capillary pistons are not autoclavable.

