

Auto-Koch

PORTABLE AUTOMATIC AUTOCLAVE





“Auto-Koch” Automatic Portable Vertical Autoclave

cylindrical chamber, single wall - die-cast alloy structure - stainless steel base - fully automatic sterilization cycle - max. operating pressure 1 atm (121 °C) - pressure gauge - automatic breather valve - automatic pressure regulator valve - safety valve - lid locking system - microprocessor controlled temperature setting system with display of set/actual temperature and sterilization time - reinforced heating element fitted on vessel bottom - chamber dimensions: Ø 310x290 h mm - overall outside dimensions Ø 440x550 h mm - power supply: 220 V, 50 Hz, 1300 W - weight 7 kg.

Code Number	Description
39159	capacity 23 litres - 1/box

The “Mediumbox” basket is used in “Auto-Koch”.

“Koch” Autoclave

semi-automatic operation - cylindrical chamber, single wall - vessel case, lid and sterilization chamber made of aluminium alloy - stainless steel base - incoloid coated heating element fitted on vessel bottom - max. working pressure 2.0 atm (121 °C) - pressure gauge - excess air breather valve - automatic pressure regulator valve - safety valve - mechanical lid locking system - inside dimensions: Ø 310x290 h mm - capacity 23 litres - overall dimensions: Ø 440x550 h mm - power supply: 220 V (110 V upon request), 50 Hz, 1300 W.



Code Number	Description
39161	capacity 23 litres - 1/box



“Field-Koch” Field Vertical Autoclave

cylindrical chamber, single wall - die-cast alloy structure - max. operating pressure 1 atm (121 °C) - pressure gauge - automatic breather valve - automatic pressure regulator valve - safety valve - lid locking system - chamber dimensions: Ø 310x295 h mm - overall outside dimensions: Ø 440x380 h mm - weight 4.5 kg.

This model is ideal for field use, where handling and independent power source become essential features. “Field-Koch” may be used with any heat source, even the naked flame. Conforms to WHO and UNICEF specifications. The “Mediumbox” basket is used in “Field-Koch”

Code Number	Description
39160	capacity 23 litres - 1/box

“Minibox” basket

For “AUTO-KOCH”, “KOCH” and “FIELD-KOCH” autoclaves - extra strong S/S perforated sheet with carrying handle

Code Number	Description
5421	dimennsions: Ø 240x210 h mm
5423	dimennsions: Ø 240x210 h mm - 4 sectors

