

Fifty 50

AUTOMATIC AUTOCLAVE CAPACITY 50 LITRES

Special rim to collect condensation when lid is opened

Manual tap for rapid steam release to accelerate cycle

Usable volume 50 litres

EC Marked

Digital microprocessor controlled panel

Certified safety valve

Pressure gauge

User friendly instruction panel according to GLP

No fitting costs: connection to the electric main is all you need

Printer (optional) for G.L.P

Drain cock for easy cleaning

Square shape for better insulation entirely manufactured in AISI 304 S/S. Easy to clean and sanitise

Electric and electronic parts comply to CEI 62-25 regulations

**10 YEARS
GUARANTEE**