

TIMO

BENCH TOP VERTICAL AUTOCLAVE CAPACITY 23 lt.



- **ABS insulated lid**
- **According to “PED” European Standard**
- **EC Marked**

Horizontal lid closure with blocking system

- *INOX Aisi 304 outer structure easy to clean and sanitize*
- *Aluminum sterilization chamber dim.: Ø 310x290 mm*
- *Overall dim.: Ø 470x580 h mm*
- *Automatic thermoregulation controlled by a microprocessor*
- *Main switch and electronic microprocessor with digital display of temperature and sterilization time*
- *Automatic air valve, electronic temperature control, electronic cut off system*
- *Reinforced heating element*
- *PT 100 accurate thermoprobe*
- *Automatic sterilization programs*
- *Maximum operating pressure: 1 ATM*
- *Maximum operating temperature: 121 °C*
- *Pressure gauge for pressure control*
- *Friendly instruction panel according to GLP and EN 17205*
- *Power supply : 220 V; 50 Hz; 1300 W; (110 V on request)*
- *Mechanical lid safety closing*
- *Portable unit*
- *IQ/OQ/PQ manuals available on request*
- *INOX basket to be ordered separately*
- *Weight: 13,250 Kg*

WHAT TO ORDER

CODE	DESCRIPTION	PACKAGING
88944	“TIMO” bench top Vertical Autoclave	1 /box
16123	“Mediumbox” basket	1 /box

Distributed By:

T 0346 26 90 90
F 0346 26 48 14
I www.salmenkipp.nl
E info@salm-en-kipp.nl