

TIMO-PRINT

BENCH TOP VERTICAL AUTOCLAVE CAPACITY 23 lt.



- ABS insulated lid
- According to "PED" European Standard
- EC Marked

Horizontal lid closure with blocking system

- INOX Aisi 304 outer structure easy to clean and sanitize
- Aluminum sterilization chamber dim.: Ø 310x290 mm
- Overall dim.: Ø 470x580 h mm
- Automatic thermoregulation controlled by a microprocessor
- Main switch and electronic microprocessor with digital display of temperature and sterilization time
- Automatic air valve, electronic temperature control, electronic cut off system
- Reinforced heating element
- PT 100 accurate thermoprobe
- Automatic sterilization programs
- Printer according to GLP
- Maximum operating pressure: 1 ATM
- Maximum operating temperature: 121 °C
- Pressure gauge for pressure control
- Friendly instruction panel according to GLP and EN 17205
- Power supply : 220 V; 50 Hz; 1300 W; (110 V on request)
- Mechanical lid safety closing
- Portable unit
- IQ/OQ/PQ manuals available on request
- INOX basket to be ordered separately
- Weight: 13,250 Kg

WHAT TO ORDER

CODE	DESCRIPTION	PACKAGING
88957	"TIMO-PRINT" bench top Vertical Autoclave	1/box
16123	"Mediumbox" basket	1/box

Distributed By:

