

Zeus

AUTOMATIC AUTOCLAVE



automatic model - cylindrical shape, single wall - entirely manufactured in AISI 304 stainless steel with polished sterilization chamber and lid - special rim placed on lid to collect condensation during lid opening - incoloid coated heating element placed on vessel bottom - PT100 precision thermoprobe made in AISI 304 stainless steel - pressure gauge - exhaust valve for fast cooling - steam safety valve - automatic valve for air dispersion - electronic cut off system - low water cut off - microprocessor controlled sterilization temperature - working pressure between 100 °C and 135 °C - max working pressure 2.1 atm - programmable sterilization time up to 99 minutes - control panel with: main switch and electronic control microprocessors with digital display of temperature and sterilization time - data recorder (optional) - capacity 140 litres - chamber dim.: 500x800 h mm - overall dimensions: 805x780x1300 h mm - weight 120 Kg - power supply: 220 V, 50 Hz, 5000 W (110V on request).

Accessories

Cat. n° **9286** - Superbox basket

Cat. n° **17494** - "Safe-lid" protection lid in s/s to be placed on the autoclave cover to protect the operator from accidental burns during sterilization cycle 1/case

Cat. n° **10171** - "Phobos" winch for easy extraction of baskets from Zeus autoclave.

Racks are not included and must be ordered separately

10
YEARS GUARANTEE

WHAT TO ORDER

CATALOGUE NUMBER	DESCRIPTION
11175	"Zeus" autoclave 1/case