

NEW

# Pipetman<sup>®</sup> Neo

## The Pipetting Standard with Lower Spring Forces

Pipetman Neo is an evolution of the Pipetman P that addresses the growing susceptibility of pipette users to repetitive strain injuries. Pipetman Neo exhibits a distinct enhancement in comfort, with the same legendary robustness, accuracy, precision and price as Pipetman P.

### Useful volume range

- Up to 50% reduction of pipetting forces.
- Same price as Pipetman P.
- Wider volume range: 10-100% of the nominal volume for the P100N and P200N models.




Easy to read volumeter

### Color-coded push buttons

- Color-coded push-buttons allow for ease of identification of the different models and corresponding quality Diamond<sup>®</sup> tips.
- Helps differentiate Pipetman Neo pipettes from the Pipetman P models.

Newly engineered springs greatly reduce susceptibility to repetitive strain injury



Color-coded push-button for easy identification

Newly engineered springs reduce susceptibility to repetitive strain injury

Pipetman Neo Pipettes				
Model	Diamond Tips	Part Number	Volume (µL)	List price
P2N 0.2 - 2 µL	D10	F144561	Min. 0.2	218,40 €
	DL10		Max. 2	
P10N 1 - 10 µL	D10	F144562	Min. 1	218,40 €
	DL10		Max. 10	
P20N 2 - 20 µL	D200	F144563	Min. 2	181,20 €
			Max. 20	
P100N 10 - 100 µL	D200	F144564	Min. 10	181,20 €
			Max. 100	
P200N 20 - 200 µL	D200	F144565	Min. 20	181,20 €
			Max. 200	
P1000N 100 - 1000 µL	D1000	F144566	Min. 100	181,20 €
			Max. 1000	



### Pipetman<sup>®</sup> Neo Starter Kits

#### Enhanced comfort. Maximum savings.

- Each Neo starter kit contains three Pipetman Neo pipettes plus boxes of certified-quality Gilson Diamond<sup>®</sup> tips. We've also included three Jimmy<sup>™</sup> microtube openers, three Single<sup>™</sup> pipette holders, a "Two Minute Inspection" poster and a Gilson Guide to Pipetting.
- Save money when compared to pipettes and accessories purchased separately.

NEW

Pipetman Neo Kits						
Model	Part Number	Pipetman Neo Models	Diamond Tips (racks of 96 tips)	Accessories	List price	Special price
Pipetman Neo Starter Kit	F167700	P20N	D200 and D1000	3 Jimmy <sup>™</sup> microtube openers 3 Single <sup>™</sup> pipette holders 1 Gilson Guide to Pipetting	525,60 €	399,00 €
		P200N				
		P100N				
Pipetman Neo PCR Starter Kit	F167600	P10N	Sterile filter tips DFL10ST, DF100ST and DF1000ST	3 Jimmy <sup>™</sup> microtube openers 3 Single <sup>™</sup> pipette holders 1 Gilson Guide to Pipetting	598,80 €	459,00 €
		P100N				
		P1000N				



Special Offer!

Purchase any Pipetman Neo Kit for a special price!  
Offer valid October 1st – December 31st 2009

Gilson maximum permissible errors are guaranteed only when Pipetman pipettes are used with the recommended Gilson Diamond<sup>®</sup> tips.

Promo Code  
KIT PROMO 09