

VARIO SHAKERS



Rack Systems for Vario Shakers


Three-Dimensional Shaker TL 10

Slow, three-dimensional orbital motion

Motion:	three-dimensional orbital motion
Shaking platform:	340 x 340 mm
Max. load:	10 kg
Shaking speed:	1 - 23 rpm
Angle of inclination:	0° to 9° (adjustable)
Timer:	0 - 120 min / continuous

Basic Equipment:

Basic device incl. basic tray and rubber mat



1 2 Rack System WS/TL
No. 0051 247

No. 6117 000

Rack Systems/Accessories:

- 1 Multi-storey rack system**
with 2 trays No. 0051 247
- 2 Additional tray for**
 - 1 size 340 x 340 mm** No. 0051 502
- 3 Additional tray for**
 - 1 size 340 x 500 mm** Nr. 0051 413

Technical Data:

Speed range:	1 - 23 rpm, steplessly adjustable
Angle of inclination:	0° - 9° (adjustable)
Max. load:	10 kg
Timer:	0 - 120 min or continuous operation
Electrical supply:	230 V or 115 V, 50/60 Hz (please indicate in case of order)
Enclosure protection:	IP 21
Heat emission	approx. 7 - 10 W
Ambient temperature:	5° C to 50° C
Relative humidity:	~ 85 %
Dimensions (wxdxh):	355 x 455 x 205 mm
Weight:	13 kg

salmenkipp

T 0346 26 90 90
F 0346 26 48 14

I www.salmenkipp.nl
E info@salm-en-kipp.nl